

 APELEM

PRÉSENTE

MOBILE RADIOGRAPHY UNITS



 **DMS**
IMAGING
A DIVISION OF THE DMS GROUP


RAFALE B
RAFALE DR
RAFALE B DR

ANALOGUE RANGE

RAFALE



RAFALE-R



RAFALE B



RAFALE

Mid-range, analogic

Power	32 kW
Voltage	125 kW
Current	400mA
Technology	Analogue/ Upgrade DR possible
Battery	No



RAFALE-R

Mid-range, analogic

Power	32 kW
Voltage	125 kV
Current	400mA
Technology	Analogue/ Upgrade DR possible
Battery	No
Rotating Arm	Yes



RAFALE B

High end, analogic, batteries

Power	32 kW / 40 kW
Voltage	125 kW
Current	320mA
Technology	Analogue/ Upgrade DR possible
Battery	Yes
Telescopic option available	Yes



UPGRADE DR

Possibility to upgrade mobile analogic in digital version with our upgrade kits – detector + tablet



iRay

Mars 1417 and 1717



Thales

Pixium EZ / DR 3543



Viewworks

Vivix 1417 and 1717



DIGITAL RANGE

RAFALE DR



RAFALE BDR



RAFALE DR

High end, digital

Power	32 kW
Voltage	125 kW
Current	400mA
Technology	Integrated DR iRay/Thales/Viewworks
Battery	No



RAFALE BDR

High end, digital, batteries

Power	32 kW / 40 kW
Voltage	125 kW
Current	320mA
Technology	Integrated DR iRay/Thales
Battery	Yes
Telescopic option available	Yes





ADDRESS

9 avenue du canal Philippe Lamour
30660 Gallargues-le-Montueux
— FRANCE

T > +33 4 67 50 49 00

F > +33 4 67 50 49 09

These products are a Class IIb medical devices, manufactured by Technix. The conformity assessment was carried out by IMQ CE0051 and bears the CE marking. It is intended for diagnosis and x-ray examinations in radiography. Read the operating instructions carefully. Photos are non-contractual.