

 APELEM

PRESENTS




MOBILE RANGE

 **DMS**
IMAGING

RAFALE
RAFALE DR

MOBILE X-RAY RANGE

ANALOGIC SOLUTION RAFALE




DIGITAL SOLUTION RAFALE DR



Mid-range, analogic

Conventional solution for customers seeking a compact, lightweight, and easy-to-operate mobile unit.

Power	32 kW
Voltage	125 kV
Current	400mA
Technology	Analogic Upgrade to digital possible
Battery	No
Rotating Arm	Option Only available with Rafale R version



UPGRADABLE TO DR

Possibility to upgrade mobile analogic in digital version with our upgrade kits – detector + tablet



iRay

Mars 1417 and 1717
DROC software



Vieworks

Vivix 1417 and 1717
VXvue Software



Mid-range, digital

Digital solution for customers seeking a compact, lightweight, and easy-to-operate mobile unit.

Power	32 kW
Voltage	125 kV
Current	400mA
Technology	Integrated Digital solution <ul style="list-style-type: none">• Pixium EZ 35x43 or 24x30• Mars 35x43
Battery	No
Rotating Arm	Option Only available with Rafale R DR version





ADDRESS > 9 Avenue du Canal Philippe LAMOUR - 30660 Gallargues-le-Montueux - FRANCE

PHONE > +33 4 66 29 09 07 – www.dms-imaging.com

These products are a Class IIb medical devices, manufactured by Technix. The conformity assessment was carried out by IMQ CE0051 and bears the CE marking. It is intended for diagnosis and x-ray examinations in radiography. Read the operating instructions carefully. Photos are non-contractual.