

**APELEM**

PRESENTS

**!M1**



## SETTING NEW STANDARDS IN MOBILE X-RAY



### OUTSTANDING MOBILITY

#### IMPROVED USABILITY

- **Smallest and lightest** fully motorized and battery operated unit
- Weight only 325 kg
- Collapsible column offers unobstructed front view

### DESIGNED FOR USERS

#### EXCEPTIONAL ERGONOMIC DESIGN

- Motorized column reduces stress and strain on technologists and radiologists
- Engineered for effortless cleaning, preventing the spread of multi-resistant bacteria

### NON-STOP OPERATION 24/7


#### ALWAYS READY FOR ACTION

- Exceptional operating time and super fast charging
- Consistently outstanding performance
- No exchange of the system battery during the life length of the unit
- Grab & go - start up time is **only 10 seconds**

## POWERED BY INTELLIGENCE

### UNCOMPROMISING EXCELLENCE IN IMAGE QUALITY

- Intuitive navigation and automated features simplifies operation allowing for quicker diagnosis
- A.I. powered image processing ensures excellent image quality, for reliable diagnosis
- Digital grid image processing enhances image rendering as well as radiographers and patients comfort



## SUSTAINABLE INVESTMENT

### RELIABLE LONG-TERM PERFORMANCE

- Outstanding expected lifespan of the X-tech cell® battery - 8-years battery warranty available
- Cost-effective solution, reducing the need for replacement and maintenance

## REVOLUTIONARY BATTERY TECHNOLOGY X-TECH CELL®

Our X-tech cell® technology keeps your mobile X-ray unit working, not waiting! We offer super fast charging and exceptional operating time. A full battery gives up to 9 hours of operating time, 1 200 exposures and 15 km in driving distance.

Competition      Operating time: **Up to 5 hours** (assuming new batteries)  
Charging from empty to full: **10 hours**



**X-tech cell®**      Operating time: **Up to 9 hours**  
Charging from empty to full: **1,5 hours**



Competition      10 minutes of charging:  
**Up to 5 minutes operating time**



**X-tech cell®**      10 minutes of charging:  
**Up to 60 minutes operating time**

