



MISTRAL

Reliable imaging, built for performance

*BRINGING INTELLIGENCE
TO IMAGING*

PRODUCT PORTFOLIO




R/F SOLUTIONS




BONE DENSITOMETRY




MOBILES



SURGICAL C-ARMS



CHEST & BONE ROOMS




RETROFIT KITS



MAMMOGRAPHY

APPLICATIONS

All types of chest and bone exams



MISTRAL: INTRODUCTION



The NEW chest & bone radiography solution

Designed for demanding environments, MISTRAL combines simplicity, durability and performance in an accessible format both for digital or analog clinical environments

MISTRAL: PRODUCT POSITIONING



- ✓ **Target users** : Clinics, hospitals, imaging centers, private practices with tight budgets
- ✓ **Applications** : Routine chest and skeletal exams
- ✓ **Our solution** : A cost-conscious and dependable manual solution to meet users' everyday needs


OPTIMIZED configuration for daily practice

- Default fixed-height table (65 cm), elevating table in option, with 4-way floating tabletop for easy positioning and patient comfort
- Column tube stand with default rotating and transverse movements
- Wide-range vertical bucky stand
- Manual and smooth movements with electromagnetic locks for safety



OPTIMIZED configuration for daily practice

- High-frequency generator up to 50 kW single or three-phase to access high image quality
- Dual-focus 300 kHU X-ray tube for both fine and large anatomical areas
- Integrated APR presets for faster workflow
- Adaptable cassette-size bucky for detectors up to 17x17"
- Optional AEC chambers and DAP



MISTRAL: MAIN ADVANTAGES

CONFIDENT cares starts with reliable system...

User-friendly interfaces and intuitive workflow

Quick and easy patient positioning

Low maintenance and robust system

Safety-focus design



ENHANCED performance with digital detectors

- Compatible with high-resolution flat panel detectors in 14x17" or 17x17" formats, wired or wireless
- All detector are managed by the same easy-to-use DR acquisition software
- Enhanced productivity with digital workflow
 - ✓ Faster image acquisition
 - ✓ Lower radiation dose with optimized image quality





THANK YOU

*BRINGING INTELLIGENCE
TO IMAGING*

 +33 4 66 29 09 07

 www.dms-imaging.com

 9 avenue du Canal Philippe Lamour
30660 Gallargues-le-Montueux
FRANCE