

APELEM

PRESENTS

!M1

ADAM  
powered by FLEXIBLE & LONG-LASTING,  
IMAGING EXCELLENCE  
EVERYWHERE**TRUE X-RAY MOBILITY®**

EFFORTLESS MOVES, EFFICIENT CARE

- **Smallest and lightest (324kg)** fully motorized for effortless transport and use
- **Reduced height** for free view forward
- **User adjustable drive handle** ensures user comfort and control

**DESIGNED FOR USERS**


COMFORT MEETS PERFORMANCE

- **Motorized column** reduces stress and strain during daily positioning
- **Smart dual workflow**, controlling key functions from both the operator screen and the collimator
- **SID guidance** at two selected distances for accurate and repeatable positioning

**POWERED BY INTELLIGENCE**

IMAGE QUALITY EXCELLENCE, SEAMLESS WORKFLOW

- **Intuitive navigation and automated features** streamline operation to support faster decisions
- **A.I. powered image processing** delivers consistent, reliable image quality for confident diagnosis
- **Digital grid image processing** enhances image rendering and improves comfort by reducing physical-grid constraints



# POWER THAT KEEPS YOU MOVING

## FAST TO CHARGE, BUILT TO LAST

Developed specifically for mobile X-ray usage, X-tech cell® technology combines high capacity with rapid charging:

- Over **1 200 exposures, 9h of use, 15km driving distance** on a single charge, for uninterrupted availability throughout the day
- A full charge takes only **1h30**, while **10 minutes of charging provides up to 1h of use** when time is critical
- **Outstanding expected lifespan** for a reliable and maintenance-free solution

## CHARGES UP TO 3X FASTER THAN COMPETITORS



OPERATION

OPERATION

OPERATION

Charging time:  
**1,5 hours**

Competition

OPERATION

CHARGING

OPERATION

Charging time:  
**5 hours**

### — !MOBILE RAD ROOM — A TURNKEY X-RAY ROOM, ANYWHERE



By combining !M1 with the **!Mobile Wall Stand**, the mobile rad room delivers the performance of a fixed X-ray room with the freedom of true mobility.



#### PLUG & PLAY

Ready to use within seconds



#### TRUE FLEXIBILITY

Easily moved between locations  
expanding imaging capacity instantly



#### COST-EFFECTIVE SOLUTION

Eliminate installation time and costs

### — !M1x —

THE GO-ANYWHERE EVOLUTION OF !M1



Building on !M1, !M1x features reinforced drive, higher clearance and vehicle-ready accessories for **fast deployment in real-world field conditions.**



#### ROBUST

Built for the outside world



#### READY TO ROLL IN SECONDS

Car loading in just 20 seconds using  
dedicated vehicle accessories



#### STAY ONLINE, RAIN OR SHINE

Connected and protected with 4G/5G  
capability and weather protection